



# ETHANOL GAS STANDARD CYLINDER REPORT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL



CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM

SFN 59282 (08/2018)

Chemical Test Operator Name (Print) <b>ANTHONY M HILL</b>		
Location <b>ROLLA</b>	Intoxilyzer® Serial Number <b>80-004204</b>	
Gas Lot Number <b>25024080A3</b>	Gas Cylinder Number <b>51</b>	Gas Expiration Date <b>10-05-2026</b>

## Check When Done:

- ☒ 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- ☒ 2. Perform an ACA Test (Level 1, Function C).  
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
  - A. 0.081 AC
  - B. 0.081 AC
  - C. 0.081 AC
- ☒ 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- ☒ 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- ☒ 5. Send the Following to the Crime Laboratory:
  - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
  - B. ACA Test Record.

Chemical Test Operator Signature 	Date <b>4-3-2025</b>
Reviewed By (Crime Laboratory Use Only) <b>Anna Narehood</b> 	Date <b>15April2025</b>

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004204  
Location = ROLL      8164.14.00 09/16  
04/03/2025      23:54

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	23:55
02 Std. Gas	0.081	23:55
03 Room Air	0.000	23:56
04 Std. Gas	0.081	23:56
05 Room Air	0.000	23:57
06 Std. Gas	0.081	23:57
07 Room Air	0.000	23:58


Lot No = 2502080A3

Cyl No = 51

Exp Date = 10/05/2026

County = 40

Oper No. = 133866

  
Operator Signature  
ANTHONY HILL

Remarks:

Form 106-I8000