



ETHANOL GAS STANDARD CYLINDER REPORT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM

SFN 59282 (08/2018)

Chemical Test Operator Name (Print) CHRISTOPHER L. ALMLIE		
Location BURK	Intoxilyzer® Serial Number 80-007097	
Gas Lot Number 25024080A3	Gas Cylinder Number 57	Gas Expiration Date 10/5/2026

Check When Done:

- ☒ 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- ☒ 2. Perform an ACA Test (Level 1, Function C).
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 - A. 0.078 AC
 - B. 0.078 AC
 - C. 0.079 AC
- ☒ 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- ☒ 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- ☒ 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature 	Date 4/7/25
Reviewed By (Crime Laboratory Use Only) Anna Narehood 	Date 16April2025

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007097
Location = BURK 8164.16.00 09/18
04/07/2025 10:57

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:59
02 Std. Gas	0.078	10:59
03 Room Air	0.000	11:00
04 Std. Gas	0.078	11:00
05 Room Air	0.000	11:00
06 Std. Gas	0.079	11:01
07 Room Air	0.000	11:01

Lot No = 25024080A3
Cyl No = 57
Exp Date = 10/05/2026
County = 07 Oper No. = 132664



Operator Signature
CHRISTOPHER ALMLIE

Remarks:

good

Form 106-I8000