



ETHANOL GAS STANDARD CYLINDER REPORT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM

SFN 59282 (08/2018)

Chemical Test Operator Name (Print) <u>GROSS, Danny J</u>		
Location <u>Lisbon PD Lib 37</u>	Intoxilyzer® Serial Number <u>80-007086</u>	
Gas Lot Number <u>25024080A3</u>	Gas Cylinder Number <u>53</u>	Gas Expiration Date <u>10/05/2026</u>

Check When Done:

- ☒ 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- ☒ 2. Perform an ACA Test (Level 1, Function C).
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
A. 0.083AC
B. 0.083AC
C. 0.083AC
- ☒ 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- ☒ 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- ☒ 5. Send the Following to the Crime Laboratory:
A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
B. ACA Test Record.

Chemical Test Operator Signature <u>Danny J Gross</u>	Date <u>04/09/2025</u>
Reviewed By (Crime Laboratory Use Only) <u>Anna Narehood</u>	Date <u>16 April 2025</u>

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007086
Location = LISB 8164.16.00 09/18
04/09/2025 18:14

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	18:15
02 Std. Gas	0.082	18:15
03 Room Air	0.000	18:16
04 Std. Gas	0.082	18:16
05 Room Air	0.000	18:17
06 Std. Gas	0.083	18:17
07 Room Air	0.000	18:18

Lot No = 25024080A3
Cyl No = 053
Exp Date = 10/05/2026
County = 37 Oper No. = 133737



Operator Signature
DANNY GROSS

Remarks:

Form 106-I8000