

**ETHANOL GAS STANDARD CYLINDER REPORT**

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM

SFN 59282 (08/2018)

Chemical Test Operator Name (Print) <b>PRESTON LANGER</b>		
Location <b>WILLISTON</b>		Intoxilyzer® Serial Number <b>80-005955</b>
Gas Lot Number <b>25024980A3</b>	Gas Cylinder Number <b>56</b>	Gas Expiration Date <b>10/5/26</b>

## Check When Done:

- ☒ 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- ☒ 2. Perform an ACA Test (Level 1, Function C).  
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).  
A. **0.081** AC  
B. **0.081** AC  
C. **0.081** AC
- ☒ 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- ☒ 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- ☒ 5. Send the Following to the Crime Laboratory:  
A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).  
B. ACA Test Record.

Chemical Test Operator Signature <b>[Signature]</b>	Date <b>4/18/25</b>
Reviewed By (Crime Laboratory Use Only) <b>Anna Narehood</b> <b>[Signature]</b>	Date <b>23April2025</b>


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005955  
Location = WILL      8164.14.00 09/16  
04/18/2025      09:22

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:23
02 Std. Gas	0.081	09:23
03 Room Air	0.000	09:24
04 Std. Gas	0.081	09:24
05 Room Air	0.000	09:25
06 Std. Gas	0.081	09:25
07 Room Air	0.000	09:25

Lot No = 25024080A3  
Cyl No = 56  
Exp Date = 10/05/2026  
County = 53      Oper No. = 130996

  
Operator Signature  
PRESTON LANGER

Remarks:

Form 106-I8000