

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80-007098 Instrument Location Tiggs 9	D					
Reason for Install/Repair						
Install After Receiving From Crime Laboratory Install After Location Change						
** Other (Specify) Install after annual inspection						
Check When Done:						
1. Surge Protector Installed/Property Grounded.						
2. Telephone Line Connected to Intoxilyzer® 8000.						
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).						
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).						
Kun Tests:						
A. Print Test (Level1, Function P).						
B-ACA Test (Level 1, Function C).						
C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio Du	ring Test).					
7. Repair and/or Maintenance Performed (if any):						
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the In	ntoxilyzer® for Use.					
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Age	ncy.					
10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair C 104-G), Print Test, ACA Test, and RFI Test.	heckout (SFN59281, Form					
Field Inspector Signature	Date					
John hele	4/23/2025					
Crime Laboratory Use Only						
This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.						
Reviewed/Certified By Anna Narehood	Certified Date 28April2025					

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007098 Location = TIOG 8164.16.00 09/18 04/23/2025 16:30

******* Printer Test ********

abcdefghijklmnopgrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA
Display Prelim Rslt? Yes Yes Display Third Digit? Inhib Printer(Y/N)? NoDisplay Volume?
Disable On Memfull? No Yes # of Print Copies? 1 Dry Select Std (D/W/I)? Standard Value? Standard Lot #? Standard Cyl #? 0.080 25024080A3 054 10/05/2026 Standard Expiration? Oper No? 130977

Flow Cal. Date: 09/26/2018 684 Slope Intercept -588446

IR Calibration Date: 09/26/2018 3um 9um

Oth Coef(*100): -40948 -15551 1st Coef(*100): 237555 2nd Coef(*100): 2571 134140 1433 H20 adj(mg/l*10k): 1097

****** Printer Test End *******

Operator Signature JOSHUA NELSON

Remarks:

I certify this is a true and accurate copy.

Form 106-I8000

Page 2 of 4

AEN

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 Location = TIOG 04/23/2025

SN 80-007098 8164.16.00 09/18 16:32

DRY CAL CHECK

Test			AC	Time
01	Room	Air	0.000	16:33
02	std.	Gas	0.081	16:33
03	Room	Air	0.000	16:34
04	std.	Gas	0.081	16:34
05	Room	Air	0.000	16:34
06	std.	Gas	0.081	16:35
07	Room	Air	0.000	16:35

Lot No = 25024080A3

Cyl No = 054

Exp Date = 10/05/2026

County = 53

Oper No. = 130977

Operator Signature JOSHUA NELSON

Remarks:

Form 106-I8000

I certify this is a true and accurate copy.

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TIOG
04/23/2025

Alcohol Analyzer SN 80-007098 8164.16.00 09/18 16:40

DRY CAL CHECK

Test			AC	Time
-	Room Room		RFI* 0.000	16:40 16:41

*Invalid Test Inhibited - RFI

Lot No = 25024080A3

Cyl No = 054

Exp Date = 10/05/2026

County = 53

Oper No. = 130977

Operator Signature JOSHUA NELSON

Remarks:

I certify this is a true and accurate copy.

Form 106-I8000