

ETHANOL GAS STANDARD CYLINDER REPORT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59282 (08/2018)

Chemical Test Operator Name (Print)	in Morten		
Location	Intoxilyzer	Intoxilyzer® Serial Number	
Gas Lot Number 08425080A1	Gas Cylinder Number	Gas Expiration Date	
Check When Done:	ו (Level 1, Function S).		
AEN \times 2. Perform an ACA Test (Level 1, Func 5/12/25 Write in the ACA Test Results (Repo A. 0.081 AC B. 0.079 AC C. 0.080 AC			
 3. Complete the Top Portion of the Into 4. File Previous Intoxilyzer® Record (S 		120-G) and Place it by the Intoxilyzer® for Use. xilvzer® Location at the Agency.	
5. Send the Following to the Crime Lab	annan o su an a straat a straat an bound an oor 🖲 saar taabaan taasaan so		
A. Ethanol Gas Standard Cylinder R B. ACA Test Record.	eport (SFN59282, Form 105-G).		
Chemical Test Operator Signature	7,4115	Date 5/7/25	

m

Reviewed By (Crime Laboratory Use Only) Anna Narehood Date

12May2025

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-004943Location = BCCC8164.14.00 09/1605/07/202516:58

	DRY CAL CHECK	
Test	AC	Time
01 Room Air	0.000	16:59
02 Std. Gas	0.081	16:59
03 Room Air	0.000	17:00
04 Std. Gas	0.079	17:00
05 Room Air	0.000	17:01
06 Std. Gas	0.080	17:01
07 Room Air	0.000	17:02

111

Operator Signature NATHAN MORTEN

Remarks:

Form 106-18000