



# ETHANOL GAS STANDARD CYLINDER REPORT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM

SFN 59282 (08/2018)

Chemical Test Operator Name (Print) <i>Brost, Easton</i>		
Location <i>DCSO</i>	Intoxilyzer® Serial Number <i>80-004947</i>	
Gas Lot Number <i>01524080A4</i>	Gas Cylinder Number <i>21</i>	Gas Expiration Date <i>3/5/26</i>

## Check When Done:

- ☒ 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- ☒ 2. Perform an ACA Test (Level 1, Function C).  
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
  - A. *0.078* AC
  - B. *0.078* AC
  - C. *0.078* AC
- ☒ 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- ☒ 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- ☒ 5. Send the Following to the Crime Laboratory:
  - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
  - B. ACA Test Record.

Chemical Test Operator Signature <i>Easton Brost</i>	Date <i>5/7/25</i>
Reviewed By (Crime Laboratory Use Only) Anna Narehood <i>Anna Narehood</i>	Date 15May2025

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004947  
Location = DUNN      8164.14.00 09/16  
05/07/2025      17:51

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	17:51
02 Std. Gas	0.078	17:52
03 Room Air	0.000	17:52
04 Std. Gas	0.078	17:53
05 Room Air	0.000	17:53
06 Std. Gas	0.078	17:53
07 Room Air	0.000	17:54

Lot No = 01524080A4

Cyl No = 21

Exp Date = 03/05/2026

County = 13

Oper No. = 132823




Operator Signature  
EASTON BROST

Remarks:

Form 106-I8000

I certify that the attached is a true and  
correct copy of the intoxilyzer test record.

Print Name/Badge #: Brost, Easton 3254

Signature: 

Date: 5/7/25 **DUNN COUNTY SHERIFF'S OFFICE**