

ETHANOL GAS STANDARD CYLINDER REPORT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59282 (08/2018)

Chemical Test Operator Name (Print) Matthew O'Brien			
Location OKPD		Intoxilyzer® Serial Number 80-007095	
Gas Lot Number 08425080A1	Gas Cylinder Number 001		Gas Expiration Date 5/5/2027

Check When Done:

- X 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).

Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).

- A. 0.079 AC
- B. 0.080 AC
- C. 0.079 AC
- X 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- X 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature			Date
That is in			5/12/2025
Reviewed By (Crime Laboratory Use Only)			Date dona coop
Anna Narehood	(1 ma	Malel	19May2025
	0-0-0-0	1	

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-007095
Location = OKPD	8164.16.00 09/18
05/12/2025	07:29

			DRY CAI	CHECK		
Т	est			AC		Time
01	Room	Air		0.000		07:30
02	Std.	Gas		0.079		07:31
03	Room	Air		0.000		07:31
04	Std.	Gas		0.080		07:31
05	Room	Air		0.000		07:32
06	Std.	Gas		0.079		07:32
07	Room	Air		0.000		07:33

Operator Signature

MATTHEW O'BRIEN

Remarks:

Form 106-18000