



**ETHANOL GAS STANDARD CYLINDER REPORT**  
NORTH DAKOTA OFFICE OF ATTORNEY GENERAL  
CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM  
SFN 59282 (08/2018)

Chemical Test Operator Name (Print) <i>Evan Wilkie</i>		
Location <i>Dickinson LEC (DKPD)</i>	Intoxilyzer® Serial Number <i>80-006689</i>	
Gas Lot Number <i>08425080A1</i>	Gas Cylinder Number <i>23</i>	Gas Expiration Date <i>05/05/2027</i>

Check When Done:

- ☒ 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- ☒ 2. Perform an ACA Test (Level 1, Function C).  
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
  - A. *0.017* AC
  - B. *0.017* AC
  - C. *0.017* AC
- ☒ 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- ☒ 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- ☒ 5. Send the Following to the Crime Laboratory:
  - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
  - B. ACA Test Record.

Chemical Test Operator Signature <i>Evan Wilkie</i>	Date <i>5/13/25</i>
Reviewed By (Crime Laboratory Use Only) Anna Narehood <i>Anna Narehood</i>	Date 20May2025

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006689  
Location = DKPD      8164.14.00 09/16  
05/13/2025      16:17

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:17
02 Std. Gas	0.077	16:18
03 Room Air	0.000	16:18
04 Std. Gas	0.077	16:18
05 Room Air	0.000	16:19
06 Std. Gas	0.077	16:19
07 Room Air	0.000	16:20

Lot No = 08425080A1  
Cyl No = 23  
Exp Date = 05/05/2027  
County = 45      Oper No. = 133098



Operator Signature  
EVAN WILKIE

Remarks:

Form 106-I8000