



ETHANOL GAS STANDARD CYLINDER REPORT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM

SFN 59282 (08/2018)

Chemical Test Operator Name (Print) <i>Tyler Rispa</i>		
Location <i>Wells CO SD</i>	Intoxilyzer® Serial Number <i>80-006491</i>	
Gas Lot Number <i>0842508041</i>	Gas Cylinder Number <i>12</i>	Gas Expiration Date <i>05/05/2027</i>

Check When Done:

- ☒ 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- ☒ 2. Perform an ACA Test (Level 1, Function C).
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 - A. *0.090* AC
 - B. *0.079* AC
 - C. *0.079* AC
- ☒ 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- ☒ 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.

AEN ☒ 5. Send the Following to the Crime Laboratory:
5/23/25 A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
B. ACA Test Record.

Chemical Test Operator Signature <i>Tyler Rispa</i>	Date <i>05/22/2025</i>
Reviewed By (Crime Laboratory Use Only) <i>Aime Nahl</i>	Date 23May2025


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006491
Location = WESO 8164.14.00 09/16
05/22/2025 09:23

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:24
02 Std. Gas	0.080	09:24
03 Room Air	0.000	09:25
04 Std. Gas	0.079	09:25
05 Room Air	0.000	09:26
06 Std. Gas	0.079	09:26
07 Room Air	0.000	09:27

Lot No = 08425080A1
Cyl No = 12
Exp Date = 05/05/2027
County = 52 Oper No. = 133682

 4771

Operator Signature
TYLER RISPA

Remarks:

Form 106-I8000