

## ETHANOL GAS STANDARD CYLINDER REPORT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59282 (08/2018)

Chemical Test Operator Name (Print) CHRISTOPHER 1,	ALMLIE	
PLPD	-	ntoxilyzer® Serial Number 80-00687
Gas Lot Number 2.5024080A3	Gas Cylinder Number 55	Gas Expiration Date

Check When Done:

1. Scan/Enter Gas Cylinder Information (Level 1, Function S).

2. Perform an ACA Test (Level 1, Function C).

Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).

- A. 0.081 AC
- В. 0.081 AC

C. 0.080 AC

X 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.

X 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.

5. Send the Following to the Crime Laboratory:

A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).

B. ACA Test Record.

Chemical Test Operator Signature	Date / ,
him	5/27/2075
Reviewed By (Crime Laboratory Use Only)	Date
Anna Narehood (Mu Mult	29May2025

FORM 105-G

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-006687Location = PLPD8164.14.00 09/1605/27/202514:07

			kinsel/				
			DRY	CAL	CHECK		
Test				AC		Time	
						_	 
01	Room	Air			0.000		14:08
02	Std.	Gas			0.081		14:08
03	Room	Air			0.000		14:09
04	Std.	Gas			0.081		14:09
05	Room	Air			0.000		14:10
06	Std.	Gas			0.080		14:10
07	Room	Air			0.000		14:11

Operator Signature

CHRISTOPHER ALMLIE

Remarks:

g 000

Form 106-18000