



ETHANOL GAS STANDARD CYLINDER REPORT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM

SFN 59282 (08/2018)

Chemical Test Operator Name (Print) <i>Chris Este</i>		
Location <i>Dickey County Courthouse</i>		Intoxilyzer® Serial Number <i>80-005453</i>
Gas Lot Number <i>08425080A1</i>	Gas Cylinder Number <i>26</i>	Gas Expiration Date <i>05-05-2022</i>

Check When Done:

- ☒ 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- ☒ 2. Perform an ACA Test (Level 1, Function C).
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
A. *0.028* AC
B. *0.028* AC
C. *0.029* AC
- ☒ 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- ☒ 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- ☒ 5. Send the Following to the Crime Laboratory:
A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
B. ACA Test Record.

Chemical Test Operator Signature <i>Chris Este</i>	Date <i>6-3-2025</i>
Reviewed By (Crime Laboratory Use Only) Anna Narehood <i>Anna Narehood</i>	Date 10June2025

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005953
Location = DCSO 8164.14.00 09/16
06/03/2025 19:04

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	19:05
02 Std. Gas	0.078	19:06
03 Room Air	0.000	19:06
04 Std. Gas	0.078	19:06
05 Room Air	0.000	19:07
06 Std. Gas	0.079	19:07
07 Room Air	0.000	19:08

Lot No = 08425080A1
Cyl No = 026
Exp Date = 05/05/2027
County = 11 Oper No. = 020504


Operator Signature
CHRIS 

Remarks:

Form 106-I8000