

**ETHANOL GAS STANDARD CYLINDER REPORT**

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM

SFN 59282 (08/2018)

Chemical Test Operator Name (Print) <i>Joshua C. McFaren</i>		
Location <i>WCPD</i>	Intoxilyzer® Serial Number <i>80-0030 60</i>	
Gas Lot Number <i>19825080A2</i>	Gas Cylinder Number <i>27</i>	Gas Expiration Date <i>9/5/2027</i>

## Check When Done:

- ☒ 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- ☒ 2. Perform an ACA Test (Level 1, Function C).  
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).  
A. *0.074* AC  
B. *0.079* AC  
C. *0.074* AC
- ☒ 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- ☒ 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- ☒ 5. Send the Following to the Crime Laboratory:  
A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).  
B. ACA Test Record.

Chemical Test Operator Signature <i>[Signature]</i>	Date <i>9/26/25</i>
Reviewed By (Crime Laboratory Use Only) <i>Anna Narehood</i>	Date <i>01Oct2025</i>

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003060  
Location = WCPD      8164.14.00 09/16  
09/26/2025      16:08

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:09
02 Std. Gas	0.079	16:09
03 Room Air	0.000	16:10
04 Std. Gas	0.079	16:10
05 Room Air	0.000	16:11
06 Std. Gas	0.079	16:11
07 Room Air	0.000	16:12


Lot No = 19825080A2

Cyl No = 27

Exp Date = 09/05/2027

County = 27

Oper No. = 131581

 #1160  
\_\_\_\_\_  
Operator Signature  
JOSHUA MCGOVEN

Remarks: "ford test"

Form 106-I8000