



ETHANOL GAS STANDARD CYLINDER REPORT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM

SFN 59282 (08/2018)

Chemical Test Operator Name (Print) <i>Allen, Kathryn</i>		
Location <i>Cass County SO-Jail Booking</i>	Intoxilyzer® Serial Number <i>80-004206</i>	
Gas Lot Number <i>01524080A4</i>	Gas Cylinder Number <i>38</i>	Gas Expiration Date <i>03-05-2026</i>

Check When Done:

- ☒ 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- ☒ 2. Perform an ACA Test (Level 1, Function C).
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 - A. 0. *083* AC
 - B. 0. *082* AC
 - C. 0. *082* AC
- ☒ 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- ☒ 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- ☒ 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature <i>Kathryn M. Allen 3671</i>	Date <i>10-22-25</i>
Reviewed By (Crime Laboratory Use Only) <i>Anna Narehood</i>	Date <i>10/23/2025</i>


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004206
Location = CASS 8164.14.00 09/16
10/22/2025 10:48

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:48
02 Std. Gas	0.083	10:49
03 Room Air	0.000	10:49
04 Std. Gas	0.082	10:49
05 Room Air	0.000	10:50
06 Std. Gas	0.082	10:50
07 Room Air	0.000	10:51

Lot No = 01524080A4
Cyl No = 38
Exp Date = 03/05/2026
County = 09 Oper No. = 090801


Operator Signature
KATHRYN ALLEN

Remarks:

All tests appear OK!

Form 106-I8000