

**ETHANOL GAS STANDARD CYLINDER REPORT**

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM

SFN 59282 (08/2018)

Chemical Test Operator Name (Print) <i>Beth Isaak</i>		
Location <i>Killdeer PD</i>		Intoxilyzer® Serial Number <i>80-004194</i>
Gas Lot Number <i>25024080A3</i>	Gas Cylinder Number <i>2</i>	Gas Expiration Date <i>10-5-26</i>

Check When Done:

☒ 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).☐ 2. Perform an ACA Test (Level 1, Function C).

Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).

A. 0. *078* ACB. 0. *078* ACC. 0. *078* AC☒ 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.☒ 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.☒ 5. Send the Following to the Crime Laboratory:

A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).

B. ACA Test Record.

Chemical Test Operator Signature <i>Beth Isaak</i>	Date <i>11-3-25</i>
Reviewed By (Crime Laboratory Use Only) <i>Anna Narehood</i>	Date <i>03Nov2025</i>

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004194
Location = KDPD 8164.14.00 09/16
11/03/2025 10:50

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:51
02 Std. Gas	0.078	10:52
03 Room Air	0.000	10:52
04 Std. Gas	0.078	10:53
05 Room Air	0.000	10:53
06 Std. Gas	0.078	10:53
07 Room Air	0.000	10:54

Lot No = 25024080A3

Cyl No = 2

Exp Date = 10/05/2026

County = 13

Oper No. = 131919



Operator Signature
BETH ISAAK

Remarks:

Form 106-I8000