



# ETHANOL GAS STANDARD CYLINDER REPORT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM

SFN 59282 (08/2018)

Chemical Test Operator Name (Print) <i>Jesse L Madche</i>		
Location <i>Nelson Co Sheriff office</i>		Intoxilyzer® Serial Number <i>80-005943</i>
Gas Lot Number <i>19825080A2</i>	Gas Cylinder Number <i>28</i>	Gas Expiration Date <i>09/05/2027</i>

## Check When Done:

- ☒ 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- ☒ 2. Perform an ACA Test (Level 1, Function C).  
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
  - A. *0.080* AC
  - B. *0.079* AC
  - C. *0.080* AC
- ☒ 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- ☒ 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- ☒ 5. Send the Following to the Crime Laboratory:
  - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
  - B. ACA Test Record.

Chemical Test Operator Signature <i>Jesse L Madche</i>	Date <i>11/5/2025</i>
Reviewed By (Crime Laboratory Use Only) <i>Anna Narehood</i>	Date <i>05Nov2025</i>

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005943  
Location = NCSO      8164.14.00 09/16  
11/05/2025      12:49

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:49
02 Std. Gas	0.080	12:50
03 Room Air	0.000	12:50
04 Std. Gas	0.079	12:51
05 Room Air	0.000	12:51
06 Std. Gas	0.080	12:52
07 Room Air	0.000	12:52

Lot No = 19825080A2  
Cyl No = 28  
Exp Date = 09/05/2027  
County = 32      Oper No. = 130917

  
Operator Signature  
JESSE MADCHE

Remarks:

*Tank Replace*  
Form 106-I8000