

ETHANOL GAS STANDARD CYLINDER REPORT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59282 (08/2018)

Chemical Test Operator Name (Print)					
pocation Intoxilyzer® Serial Number 80-0000					
Gas Lot Number Gas Cylinder	r Number 018	Gas Expiration Date 09/5/2027			
Check When Done:		,			
1. Scan/Enter Gas Cylinder Information (Level 1, Fur	nction S).				
2. Perform an ACA Test (Level 1, Function C).					
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).					
А. 0 <i>Q</i> \$ ♥ АС В. 0.0 В Q АС С. 0. 0 8 Q АС					
3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.					
💢 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.					
5. Send the Following to the Crime Laboratory:					
A. Ethanol Gas Standard Cylinder Report (SFN59 B. ACA Test Record.	9282, Form 105-G).				
Chemiçal Test Operator Signature		Date			
Criselda A. Crum		11-21-25			
Reviewed By (Grime Laboratory Use Only)	n Anna Inge	emanson 11-28-25			

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = MOSO
11/21/2025

Alcohol Analyzer SN 80-006676 8164.14.00 09/16 22:27

DRY CAL CHECK

Test			AC	Time
01	Room	Air	0.000	22:28
02	Std.	Gas	0.080	22:28
03	Room	Air	0.000	22:29
04	Std.	Gas	0.080	22:29
05	Room	Air	0.000	22:30
06	Std.	Gas	0.080	22:30
07	Room	Air	0.000	22:30

Lot No = 19825080A2

Cyl No = 018

Exp Date = 09/05/2027

County = 31

Oper No. = 130485

Operator Signature CRISELDA CRUZ

Remarks:

Form 106-I8000

The undersigned, having custody of the original public record, certifies that the information hereon is a true and correct copy of the original document maintained as part of the activities of this office.