



# ETHANOL GAS STANDARD CYLINDER REPORT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM

SFN 59282 (08/2018)

Chemical Test Operator Name (Print) <b>CHRISTOPHER ALMLIE</b>		
Location <b>PLPD</b>		Intoxilyzer® Serial Number <b>80-006687</b>
Gas Lot Number <b>25024080A3</b>	Gas Cylinder Number <b>55</b>	Gas Expiration Date <b>10/05/2026</b>

## Check When Done:

- ☒ 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- ☒ 2. Perform an ACA Test (Level 1, Function C).  
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
  - A. **0.082** AC
  - B. **0.081** AC
  - C. **0.081** AC
- ☒ 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- ☒ 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- ☒ 5. Send the Following to the Crime Laboratory:
  - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
  - B. ACA Test Record.

Chemical Test Operator Signature 	Date <b>11/25/25</b>
Reviewed By (Crime Laboratory Use Only) <b>Anna Ingemansen</b> 	Date <b>06Jan2026</b>

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006687  
Location = PLPD      8164.14.00 09/16  
11/25/2025      14:46

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:47
02 Std. Gas	0.082	14:48
03 Room Air	0.000	14:48
04 Std. Gas	0.081	14:49
05 Room Air	0.000	14:49
06 Std. Gas	0.081	14:49
07 Room Air	0.000	14:50

Lot No = 25024080A3  
Cyl No = 55  
Exp Date = 10/05/2026  
County = 07      Oper No. = 132664



Operator Signature  
CHRISTOPHER ALMLIE

Remarks:

good

Form 106-I8000