



**INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT**  
 NORTH DAKOTA OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM  
 SFN 59281 (06/2018)

Serial Number	80-006679	Instrument Location	HILL
Reason for Install/Repair			
<input checked="" type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input checked="" type="checkbox"/> Other (Specify) <u>Install after breath tube repair</u>			

Check When Done:

- 1. Surge Protector Installed/Property Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
  - A. Print Test (Level 1, Function P).
  - B. ACA Test (Level 1, Function C).
  - C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
- 7. Repair and/or Maintenance Performed (if any): NA
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature	<i>CWB</i>	Date	2-12-2026
---------------------------	------------	------	-----------

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Reviewed/Certified By	Anna Ingemansen	<i>Anna Ingemansen</i>	Certified Date	17Feb2026
-----------------------	-----------------	------------------------	----------------	-----------



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006679  
Location = HILL      8164.14.00 09/16  
02/12/2026      17:27

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	17:28
02 Std. Gas	0.080	17:28
03 Room Air	0.000	17:29
04 Std. Gas	0.081	17:29
05 Room Air	0.000	17:30
06 Std. Gas	0.080	17:30
07 Room Air	0.000	17:31

Lot No = 08425080A1  
Cyl No = 15  
Exp Date = 05/05/2027  
County = 49      Oper No. = 131630

  
\_\_\_\_\_  
Operator Signature  
AUSTIN BARIL

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006679  
Location = HILL                    8164.14.00 09/16  
02/12/2026                            17:31

Test	AC	Time
01 Room Air	0.000	17:33
02 *Subject Test	RFI*	17:34
03 Room Air	0.000	17:34

\*Invalid Test  
  Inhibited - RFI

Sub Name = TEST, SUBJECT NA  
Sub DOB = 01/01/2000  
Sub Sex = Male                      Weight = NA  
Test = OTH                            Cit = NA  
Dr. Lic. = ND/NA  
Lot No = 08425080A1  
Cyl No = 15  
Expiration Date = 05/05/2027  
County = 49                          Oper No. = 131630



Operator Signature  
AUSTIN BARIL

Remarks:

Form 106-I8000