



ETHANOL GAS STANDARD CYLINDER REPORT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL


CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM

SFN 59282 (08/2018)

Chemical Test Operator Name (Print) <i>PRESTON LANGER</i>		
Location <i>WILLESTON</i>	Intoxilyzer® Serial Number <i>80-005939</i>	
Gas Lot Number <i>19825060A2</i>	Gas Cylinder Number <i>30</i>	Gas Expiration Date <i>9/5/27</i>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 - A. *0.080* AC
 - B. *0.081* AC
 - C. *0.080* AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature 	Date <i>3/2/26</i>
Reviewed By (Crime Laboratory Use Only) Anna Ingemansen 	Date 04March2026


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005939
Location = WILL 8164.14.00 09/16
03/02/2026 10:12

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:13
02 Std. Gas	0.080	10:13
03 Room Air	0.000	10:14
04 Std. Gas	0.081	10:14
05 Room Air	0.000	10:14
06 Std. Gas	0.080	10:15
07 Room Air	0.000	10:15

Lot No = 19825080A2
Cyl No = 30
Exp Date = 09/05/2027
County = 53 Oper No. = 130996



Operator Signature
PRESTON LANGER

Remarks:

Form 106-I8000