



ETHANOL GAS STANDARD CYLINDER REPORT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM

SFN 59282 (08/2018)

Chemical Test Operator Name (Print) <i>Daniel Magelky</i>		
Location <i>BMDC</i>	Intoxilyzer® Serial Number <i>40-004936</i>	
Gas Lot Number <i>26-4437</i>	Gas Cylinder Number <i>113</i>	Gas Expiration Date <i>1/29/2028</i>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 - A. *0078* AC
 - B. *0078* AC
 - C. *0078* AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature <i>Daniel Magelky</i>	Date <i>4/21/26</i>
Reviewed By (Crime Laboratory Use Only) <i>Anna Ingemansen</i> Anna Ingemansen	Date 30April2026

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004936
Location = BMDC 8164.14.00 09/16
04/21/2026 09:48

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:49
02 Std. Gas	0.078	09:49
03 Room Air	0.000	09:50
04 Std. Gas	0.078	09:50
05 Room Air	0.000	09:51
06 Std. Gas	0.078	09:51
07 Room Air	0.000	09:52

Lot No = 264437
Cyl No = 113
Exp Date = 01/29/2028
County = 08 Oper No. = 133598



Operator Signature
DANIEL MAGELKY

Remarks:

Form 106-I8000