



# ETHANOL GAS STANDARD CYLINDER REPORT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM

SFN 59282 (08/2018)

|   |  |  |
|---|--|--|
| Chemical Test Operator Name (Print)<br><i>Alex Cook</i> |  |  |
| Location<br><i>MOHL</i>                                 | Intoxilyzer® Serial Number<br><i>80-006689</i> |  |
| Gas Lot Number<br><i>08425080A1</i>                     | Gas Cylinder Number<br><i>013</i>              | Gas Expiration Date<br><i>05/05/27</i> |

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).  
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
  - A. 0.*073* AC
  - B. 0.*078* AC
  - C. 0.*078* AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
  - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
  - B. ACA Test Record.

|   |                        |
|---|------------------------|
| Chemical Test Operator Signature<br><i>[Signature]</i>                            | Date<br><i>6/10/26</i> |
| Reviewed By (Crime Laboratory Use Only)<br>Anna Ingemansen <i>Anna Ingemansen</i> | Date<br>11 June 2026   |

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006689  
Location = MOHL      8164.14.00 09/16  
06/10/2026      08:46

DRY CAL CHECK

| Test        | AC    | Time  |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 08:47 |
| 02 Std. Gas | 0.078 | 08:47 |
| 03 Room Air | 0.000 | 08:48 |
| 04 Std. Gas | 0.078 | 08:48 |
| 05 Room Air | 0.000 | 08:49 |
| 06 Std. Gas | 0.078 | 08:49 |
| 07 Room Air | 0.000 | 08:50 |

Lot No = 08425080A1  
Cyl No = 013  
Exp Date = 05/05/2027  
County = 38      Oper No. = 133141

  
\_\_\_\_\_  
Operator Signature  
ALEX COOK

Remarks:

Form 106-I8000