



ETHANOL GAS STANDARD CYLINDER REPORT
NORTH DAKOTA OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM
SFN 59282 (08/2018)

Chemical Test Operator Name (Print) <i>Hauf, Anthony, E</i>		
Location <i>Rugby</i>	Intoxilyzer® Serial Number <i>80-004953</i>	
Gas Lot Number <i>264944</i>	Gas Cylinder Number <i>312</i>	Gas Expiration Date <i>5/12/2028</i>

Check When Done:

- ☒ 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- ☒ 2. Perform an ACA Test (Level 1, Function C).
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 - A. *0.081* AC
 - B. *0.081* AC
 - C. *0.080* AC
- ☒ 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- ☒ 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- ☒ 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature <i>Anthony Hauf</i>	Date <i>9/9/26</i>
Reviewed By (Crime Laboratory Use Only) <i>Anna Ingemansen</i>	Date 10Sept2026

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004953
Location = RUGB 8164.14.00 09/16
09/09/2026 16:46

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:48
02 Std. Gas	0.081	16:48
03 Room Air	0.000	16:49
04 Std. Gas	0.081	16:49
05 Room Air	0.000	16:50
06 Std. Gas	0.080	16:50
07 Room Air	0.000	16:51

Lot No = 264944

Cyl No = 312

Exp Date = 05/12/2028

County = 35

Oper No. = 132732



#4525

Operator Signature
ANTHONY HAUF

Remarks:

Form 106-I8000