



# ETHANOL GAS STANDARD CYLINDER REPORT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM

SFN 59282 (08/2018)

|                                                          |                                             |                                      |
|----------------------------------------------------------|---------------------------------------------|--------------------------------------|
| Chemical Test Operator Name (Print) <u>Luke Olmstead</u> |                                             |                                      |
| Location <u>Pembina County Jail</u>                      | Intoxilyzer® Serial Number <u>80-006494</u> |                                      |
| Gas Lot Number <u>264944</u>                             | Gas Cylinder Number <u>301</u>              | Gas Expiration Date <u>5/12/2028</u> |

## Check When Done:

- ☒ 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- ☒ 2. Perform an ACA Test (Level 1, Function C).  
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
  - A. 0.080 AC
  - B. 0.080 AC
  - C. 0.080 AC
- ☒ 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- ☒ 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- ☒ 5. Send the Following to the Crime Laboratory:
  - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
  - B. ACA Test Record.

|                                                                                       |                        |
|---------------------------------------------------------------------------------------|------------------------|
| Chemical Test Operator Signature <u>Luke Olmstead</u>                                 | Date <u>9/18/2026</u>  |
| Reviewed By (Crime Laboratory Use Only) <u>Anna Ingemansen</u> <u>Anna Ingemansen</u> | Date <u>21Sept2026</u> |

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006494  
Location = CAVA      8164.14.00 09/16  
09/18/2026      09:29

DRY CAL CHECK

| Test        | AC    | Time  |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 09:30 |
| 02 Std. Gas | 0.080 | 09:30 |
| 03 Room Air | 0.000 | 09:31 |
| 04 Std. Gas | 0.080 | 09:31 |
| 05 Room Air | 0.000 | 09:32 |
| 06 Std. Gas | 0.080 | 09:32 |
| 07 Room Air | 0.000 | 09:33 |

Lot No = 264944  
Cyl No = 301  
Exp Date = 05/12/2028  
County = 34      Oper No. = 134363



Operator Signature  
LUKE OLMSTEAD

Remarks:

Form 106-I8000