



# ETHANOL GAS STANDARD CYLINDER REPORT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM

SFN 59282 (08/2018)

|                                                              |                                                |                                       |
|--------------------------------------------------------------|------------------------------------------------|---------------------------------------|
| Chemical Test Operator Name (Print)<br><i>Tara Cavanaugh</i> |                                                |                                       |
| Location<br><i>BIA-Spirit Lake Police</i>                    | Intoxilyzer® Serial Number<br><i>80-085362</i> |                                       |
| Gas Lot Number<br><i>26-4437</i>                             | Gas Cylinder Number<br><i>#110</i>             | Gas Expiration Date<br><i>1-29-28</i> |

## Check When Done:

- ☒ 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- ☒ 2. Perform an ACA Test (Level 1, Function C).  
Write in the ACA Test Results (Reported to 3 Digits; Ex. 0.081).
  - A. *0.078* AC
  - B. *0.079* AC
  - C. *0.079* AC
- ☒ 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- ☒ 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- ☒ 5. Send the Following to the Crime Laboratory:

AEI  
25Sep26

- A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
- B. ACA Test Record.

|                                                                   |                          |
|-------------------------------------------------------------------|--------------------------|
| Chemical Test Operator Signature<br><i>Tara Cavanaugh</i>         | Date<br><i>9-10-26</i>   |
| Reviewed By (Crime Laboratory Use Only)<br><i>Anna Ingemansen</i> | Date<br><i>25Sep2026</i> |

Checklist #6 was checked by Anna Ingemansen. The two tests were received at the Crime Laboratory. This casues no issues with the Ethanol Gas Standard Cylinder Report and Cylinder Installation.  
AEI 25Sep2026

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005362  
Location = FTOT      8164.14.00 09/16  
09/06/2026      13:04

DRY CAL CHECK

| Test        | AC           | Time  |
|-------------|--------------|-------|
| 01 Room Air | 0.000        | 13:05 |
| 02 Std. Gas | <u>0.078</u> | 13:05 |
| 03 Room Air | 0.000        | 13:06 |
| 04 Std. Gas | <u>0.079</u> | 13:06 |
| 05 Room Air | 0.000        | 13:07 |
| 06 Std. Gas | <u>0.079</u> | 13:07 |
| 07 Room Air | 0.000        | 13:08 |


Lot No = 264437

Cyl No = 110

Exp Date = 01/29/2028

County = 62

Oper No. = 133407



Operator Signature  
TARA CAVANAUGH

Remarks:

Form 106-I8000