



ETHANOL GAS STANDARD CYLINDER REPORT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM

SFN 59282 (05/2026)

Chemical Test Operator Name (Print) <i>Nickolas Oster</i>		
Location <i>Tioga PD (TI06)</i>		Intoxilyzer® Serial Number <i>80-007098</i>
Gas Lot Number <i>264944</i>	Gas Cylinder Number <i>295</i>	Gas Expiration Date <i>5/12/2028</i>

Check When Done:

- ☒ 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- ☒ 2. Perform an ACA Test (Level 1, Function C).
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
A. 0. 0 8 0 AC
B. 0. 0 8 1 AC
C. 0. 0 8 1 AC
- ☒ 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- ☒ 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- ☒ 5. Send the Following to the Crime Laboratory:
A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
B. ACA Test Record.

Chemical Test Operator Signature <i>Nickolas Oster</i>	Date <i>10/2/20</i>
---	------------------------

Crime Laboratory Use Only

Reviewed By (Print Name) Anna Ingemansen	Reviewed By Signature <i>Anna Ingemansen</i>	Date 05Oct2026
---	---	-------------------

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007098
Location = TIOG 8164.16.00 09/18
10/02/2026 18:46

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	18:47
02 Std. Gas	0.080	18:47
03 Room Air	0.000	18:48
04 Std. Gas	0.081	18:48
05 Room Air	0.000	18:49
06 Std. Gas	0.081	18:49
07 Room Air	0.000	18:50

Lot No = 264944
Cyl No = 295
Exp Date = 05/12/2028
County = 53 Oper No. = 131838



Operator Signature
NICKOLAS OSTER

Remarks:

ACA Test

Form 106-I8000