



ETHANOL GAS STANDARD CYLINDER REPORT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM

SFN 59282 (08/2018)

Chemical Test Operator Name (Print) <i>Terry Dillon</i>		
Location <i>CARS</i>	Intoxilyzer® Serial Number <i>80-006681</i>	
Gas Lot Number <i>264944</i>	Gas Cylinder Number <i>303</i>	Gas Expiration Date <i>5/12/28</i>

Check When Done:

- ☒ 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- AEI ☒ 2. Perform an ACA Test (Level 1, Function C).
10/5/26 Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
A. 0.079 AC
B. 0.079 AC
C. 0.079 AC
- ☒ 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- ☒ 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- ☒ 5. Send the Following to the Crime Laboratory:
A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
B. ACA Test Record.

Chemical Test Operator Signature <i>[Signature]</i>	Date <i>10/3/26</i>
Reviewed By (Crime Laboratory Use Only) Anna Ingemansen <i>Anna Ingemansen</i>	Date 05Oct2026

Anna Ingemansen checked off #2 as the ACA test is attached with the test results written.

AEI 10/5/2026

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006681
Location = CARS 8164.14.00 09/16
10/03/2026 18:08

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	18:09
02 Std. Gas	0.078	18:09
03 Room Air	0.000	18:10
04 Std. Gas	0.078	18:10
05 Room Air	0.000	18:10
06 Std. Gas	0.079	18:11
07 Room Air	0.000	18:11

Lot No = 264944
Cyl No = 303
Exp Date = 05/12/2028
County = 19 Oper No. = 132318



Operator Signature
TERRY DILLON

Remarks:

Form 106-I8000